

CD81 Hamster anti-Mouse Monoclonal (FITC) Antibody - LS-C5608 - LSBio

CatalogID:	LS-C5608
Target:	CD81 molecule
Synonyms:	CD81 Antibody, CD81 antigen Antibody, CD81 molecule Antibody, CVID6 Antibody, S5.7 Antibody, TAPA1 Antibody, TAPA-1 Antibody, Tetraspanin-28 Antibody, TSPAN28 Antibody, Tspan-28 Antibody
Family / Subfamily:	Tetraspan / not assigned-Tetraspan
Host	CD81 antibody was produced in Hamster
Clonality:	Monoclonal
Isotype:	IgG3,k
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD81 antibody was raised against Mouse
Immunogen:	CD81 antibody was raised against mouse epithelial cell line PAM212.
Specificity:	Recognizes mouse CD81 (Mr 25kD).
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, sodium azide.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow Cytometry: 1 ug/1E6 cells.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences