

Listeria monocytogenes Rabbit anti-Listeria monocytogenes Polyclonal (FITC) Antibody - LS-C56360 - LSBio

CatalogID:	LS-C56360
Target:	Listeria monocytogenes
Host	Rabbit
Clonality:	Polyclonal
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Listeria monocytogenes
Antigen Type:	Cells
Immunogen:	Listeria monocytogenes; ATCC strain 43251.
Specificity:	whole cells
Reactivity:	Listeria monocytogenes
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, 10mg/ml BSA
Recommended Storage:	Store at +4C or -20C if preferred. Product should be stored undiluted. Storage in frost free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freeze/thaw cycles.
Usage Summary:	Immunohistology: Acetone is recommended.
Uses:	Immunofluorescence (1:10 - 1:50) (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Concentration:	4 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences