

**CD86 Rat anti-Mouse Monoclonal (PE) Antibody - LS-C5638 - LSBio**

<b>CatalogID:</b>	LS-C5638
<b>Target:</b>	CD86 molecule
<b>Synonyms:</b>	CD86 Antibody, Activation B7-2 antigen Antibody, B7-2 Antibody, B7.2 Antibody, B70 Antibody, CD86 molecule Antibody, CTLA-4 counter-receptor B7.2 Antibody, CD86 antigen Antibody, FUN-1 Antibody, LAB72 Antibody, BU63 Antibody, CD28LG2 Antibody, Ctla-4 counter-receptors b7.2 Antibody
<b>Host</b>	CD86 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2b,k
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	CD86 antibody was raised against Mouse
<b>Specificity:</b>	Recognizes mouse CD86.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 0.1% sodium azide, stabilizing agent.
<b>Recommended Storage:</b>	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
<b>Usage Summary:</b>	Identification and enumeration of CD86+ cells by flow cytometry. Identification of CD86+ cells in frozen sections. Immunoprecipitation. In vivo and in vitro blocking assays. Flow Cytometry: < 0.2 ug/1E6 cells.
<b>Uses:</b>	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences