

CD86 Rat anti-Mouse Monoclonal (PE, Cy5) Antibody - LS-C5666 - LSBio

CatalogID:	LS-C5666
Target:	CD86 molecule
Synonyms:	CD86 Antibody, Activation B7-2 antigen Antibody, B7-2 Antibody, B7.2 Antibody, B70 Antibody, CD86 molecule Antibody, CTLA-4 counter-receptor B7.2 Antibody, CD86 antigen Antibody, FUN-1 Antibody, LAB72 Antibody, BU63 Antibody, CD28LG2 Antibody, Ctla-4 counter-receptors b7.2 Antibody
Host	CD86 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Conjugations:	Phycoerythrin (PE), Cy5
Immunogen Species:	CD86 antibody was raised against Mouse
Immunogen:	CD86 antibody was raised against IPS-activated mouse B cells.
Specificity:	Recognizes mouse CD86, the B7-2 co-stimulatory molecule and a ligand for CD28 and CD152/CTLA-4.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. In vivo and in vitro blocking assays. Flow Cytometry: 0.1 ug/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences