

THY1 / CD90 Rat anti-Mouse Monoclonal Antibody - LS-C5683 - LSBio	
CatalogID:	LS-C5683
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG1
Immunogen Species:	THY1 / CD90 antibody was raised against Mouse
Antigen Type:	Cells
Immunogen:	THY1 / CD90 antibody was raised against purified mouse LN T cells.
Specificity:	Recognizes mouse CD90 Thy-1.
Reactivity:	Mouse
Purification:	Ion exchange chromatography
Reconstitution:	deionized water. Possible additional volumes for resuspension: 500 µl
Presentation:	Lyophilized.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Cytotoxic and staining antibody. For staining, use 10 ul/1E6 cells. Exact titer for cytotoxicity should be determined by the investigator.
Uses:	Cytotoxic Assay (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	