

THY1 / CD90 Mouse anti-Human Monoclonal (FITC) Antibody - LS-C5695 - LSBio

CatalogID:	LS-C5695
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	THY1 / CD90 antibody was raised against Human
Antigen Type:	Purified protein
Immunogen:	THY1 / CD90 antibody was raised against purified human brain Thy-1.
Specificity:	Recognizes the human CD90 cell surface antigen. Species cross-reactivity: Cynomolgus monkey.
Reactivity:	Human, Monkey
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 1% BSA, 0.1% sodium azide.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Flow Cytometry. Flow Cytometry: Neat. Use 10 ul to label 1E6 cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences