

THY1 / CD90 Mouse anti-Rat Monoclonal (PE) Antibody - LS-C5702 - LSBio

CatalogID:	LS-C5702
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	THY1 / CD90 antibody was raised against Rat
Immunogen:	THY1 / CD90 antibody was raised against rat Thy-1.
Specificity:	Recognizes rat and mouse CD90, also known as Thy1.1, expressed on a variety of cell types including thymocytes, neuronal cells, stem cells, T lymphocytes, immature B cells and connective tissues.
Reactivity:	Rat, Mouse
Purification:	Ion exchange chromatography
Reconstitution:	deionized water
Presentation:	Lyophilized.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Flow Cytometry. Flow Cytometry: Neat. Use 10 ul to label 1E6 cells in 100ul.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences