

THY1 / CD90 Rat anti-Mouse Monoclonal Antibody - LS-C5714 - LSBio	
<b>CatalogID:</b>	LS-C5714
<b>Target:</b>	Thy-1 cell surface antigen (THY1)
<b>Synonyms:</b>	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
<b>Host</b>	THY1 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Immunogen Species:</b>	THY1 / CD90 antibody was raised against Mouse
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	THY1 / CD90 antibody was raised against purified mouse LN T cells.
<b>Specificity:</b>	Recognizes mouse CD90 Thy-1.1.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Ion exchange chromatography
<b>Reconstitution:</b>	deionized water. Possible additional volumes for resuspension: 500 µl
<b>Presentation:</b>	Lyophilized.
<b>Recommended Storage:</b>	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
<b>Usage Summary:</b>	Suitable for use in Immunocytochemistry as a cytotoxic and staining antibody. Extremely efficient for killing in vivo. Immunocytochemistry: For staining, use 10 ul/1E6 cells.
<b>Uses:</b>	ICC (Optimal dilution to be determined by the researcher)
<b>Size:</b>	500 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	