

THY1 / CD90 Rat anti-Mouse Monoclonal (G7) Antibody - LS-C5720 - LSBio

CatalogID:	LS-C5720
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2c,k
Clone Name:	G7
Immunogen Species:	THY1 / CD90 antibody was raised against Mouse
Specificity:	Recognizes Thy-1.1 and Thy-1.2 alloantigens (Thy-1 epitope region A).
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS. No preservative added.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunocytochemistry (acetone-fixed, frozen sections). In vitro T cell activation. Induction of apoptosis.
Uses:	ICC, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences