

THY1 / CD90 Rat anti-Mouse Monoclonal (PE, Cy5) (G7) Antibody - LS-C5725 - LSBio

CatalogID:	LS-C5725
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2c,k
Clone Name:	G7
Conjugations:	Phycoerythrin (PE), Cy5
Immunogen Species:	THY1 / CD90 antibody was raised against Mouse
Specificity:	Recognizes Thy-1.1 and Thy-1.2 alloantigens (Thy-1 epitope region A).
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Do not freeze. Stable at 4° as an undiluted liquid. Dilute only prior to immediate use. Freezing PE conjugates will result in a substantial loss of enzymatic activity.
Usage Summary:	Suitable for use in Flow Cytometry, in vitro T cell activation and Induction of apoptosis. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	0.1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences