

**THY1 / CD90 Rat anti-Mouse Monoclonal (APC) (30-H12) Antibody - LS-C5731 - LSBio**

<b>CatalogID:</b>	LS-C5731
<b>Target:</b>	Thy-1 cell surface antigen (THY1)
<b>Synonyms:</b>	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
<b>Host</b>	THY1 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2b,k
<b>Clone Name:</b>	30-H12
<b>Conjugations:</b>	Allophycocyanin (APC)
<b>Immunogen Species:</b>	THY1 / CD90 antibody was raised against Mouse
<b>Specificity:</b>	Recognizes mouse CD90.2/Thy-1.2 alloantigen, Mr 18kD.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, sodium azide and a stabilizing agent
<b>Recommended Storage:</b>	Store at 4°C. Do not freeze.
<b>Usage Summary:</b>	Suitable for use in Flow Cytometry, Immunohistochemistry and. In vivo and in vitro depletion of T cells. Flow Cytometry: 0.2 ug/1E6 cells. Immunohistochemistry: Frozen.
<b>Uses:</b>	IHC - Frozen, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Concentration:</b>	0.1 mg/ml
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences