

THY1 / CD90 Rat anti-Mouse Monoclonal (PE, Cy5) (30-H12) Antibody - LS-C5732 - LSBio

CatalogID:	LS-C5732
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2b,k
Clone Name:	30-H12
Conjugations:	Phycoerythrin (PE), Cy5
Immunogen Species:	THY1 / CD90 antibody was raised against Mouse
Specificity:	Recognizes mouse CD90.2/Thy-1.2 alloantigen, Mr 18kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunohistochemistry (frozen sections). In vivo and in vitro depletion of T cells. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	IHC - Frozen, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences