

**ICAM2 / CD102 Rat anti-Mouse Monoclonal (PE) Antibody - LS-C5763 - LSBio**

<b>CatalogID:</b>	LS-C5763
<b>Target:</b>	intercellular adhesion molecule 2 (ICAM2)
<b>Synonyms:</b>	ICAM2 Antibody, CD102 Antibody, ICAM-2 Antibody, CD102 antigen Antibody
<b>Family / Subfamily:</b>	ICAM / not assigned-ICAM
<b>Host</b>	ICAM2 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a,k
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	ICAM2 / CD102 antibody was raised against Mouse
<b>Immunogen:</b>	ICAM2 / CD102 antibody was raised against cOS cells transfected with mouse ICAM-2 cDNA.
<b>Specificity:</b>	Reacts with mouse ICAM-2 (CD102) cell surface glycoprotein.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 0.1% sodium azide, stabilizing agent.
<b>Recommended Storage:</b>	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
<b>Usage Summary:</b>	Flow cytometry. In vitro blocking assays. Immunohistochemistry (frozen sections). Immunoprecipitation. Flow Cytometry: 0.2 ug/1E6 cells.
<b>Uses:</b>	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences