

**VCAM1 / CD106 Mouse anti-Human Monoclonal (PE) Antibody - LS-C5782 - LSBio**

<b>CatalogID:</b>	LS-C5782
<b>Target:</b>	vascular cell adhesion molecule 1 (VCAM1)
<b>Synonyms:</b>	VCAM1 Antibody, CD106 Antibody, CD106 antigen Antibody, INCAM-100 Antibody, V-CAM 1 Antibody, VCAM-1 Antibody
<b>Family / Subfamily:</b>	ICAM / not assigned-ICAM
<b>Host</b>	VCAM1 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	VCAM1 / CD106 antibody was raised against Human
<b>Specificity:</b>	Recognizes human CD106/VCAM-1.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 0.1% sodium azide, stabilizing agent.
<b>Recommended Storage:</b>	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
<b>Usage Summary:</b>	Suitable for use in Flow Cytometry, Immunohistochemistry, ELISA and Western Blot. Flow Cytometry: 10 ul/1E6 cells. Immunohistochemistry: frozen sections.
<b>Uses:</b>	IHC - Frozen, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences