

MHC Class II RT-1B Mouse anti-Rat Monoclonal (FITC) (OX-6) Antibody - LS-C58496 - LSBio

CatalogID:	LS-C58496
Target:	MHC Class II RT-1B
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	OX-6
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Rat
Antigen Type:	Membrane extract
Immunogen:	Rat thymocyte membrane glycoproteins.
Specificity:	Recognizes a monomorphic determinant of the rat RT1B MHC class II antigen present on B lymphocytes, dendritic cells, some macrophages and certain epithelial cells. This antibody Reacts with certain mouse strains of MHC haplotypes k and s. The clone does not react with the Rat BDIX strain.
Reactivity:	Rat, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 1% BSA.
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Usage Summary:	Flow Cytometry: Use 10 ul of the suggested working dilution to label 10 ⁶ cells in 100 ul. Method sheets are available on request.
Uses:	Flow Cytometry (1:1 - 1:10) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences