

Glucose Oxidase Sheep anti-Aspergillus niger Polyclonal (TR) Antibody - LS-C59327 - LSBio	
CatalogID:	LS-C59327
Target:	Glucose Oxidase
Host	Glucose Oxidase antibody was produced in Sheep
Clonality:	Polyclonal
Conjugations:	Texas Red (TR)
Immunogen Species:	Glucose Oxidase antibody was raised against Aspergillus niger
Antigen Type:	Purified protein
Immunogen:	Glucose Oxidase antibody was raised against glucose Oxidase [Aspergillus niger].
Specificity:	Glucose Oxidase [Aspergillus niger]. Cross reactivity against Glucose Oxidase from other tissues and species may occur but have not been specifically determined
Reactivity:	Aspergillus niger
Purification:	Purified IgG
Presentation:	0.01 M sodium phosphate, 0.14 M sodium chloride, pH 7.4.
Recommended Storage:	Store at 4°C.
Usage Summary:	Suitable for immunoblotting (western or dot blot), ELISA, immunoprecipitation and most immunological methods requiring high titer and specificity. This product has been assayed against 1.0 ug of Glucose Oxidase [Aspergillus niger] in a standard sandwich ELISA using Peroxidase conjugated Affinity Purified anti- Sheep IgG [H&L] (Sheep) catalog no. LS-C60884 and ABTS as a substrate for 30 minutes at room temperature.
Uses:	Western blot, Immunoprecipitation, ELISA (1:1000 - 1:5000) (Optimal dilution to be determined by the researcher)
Size:	2,000 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	