

Syrian Hamster IgG Rabbit anti-Syrian hamster Polyclonal (TR) Antibody - LS-C59363 - LSBio

CatalogID:	LS-C59363
Target:	Syrian Hamster IgG
Host	Rabbit
Clonality:	Polyclonal
Conjugations:	Texas Red (TR)
Immunogen Species:	Syrian hamster
Antigen Type:	Purified protein
Immunogen:	Hamster IgG whole molecule.
Specificity:	Hamster IgG and Hamster Serum
Reactivity:	Syrian hamster, Hamster
Purification:	Purified IgG
Presentation:	0.01 M sodium phosphate, 0.14 M sodium chloride, pH 7.4.
Recommended Storage:	Store at 4°C.
Usage Summary:	Suitable for classical and molecular immunological methods.
Uses:	Undetermined
Size:	2,000 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences