

Human IgG Fab'2 Rabbit anti-Human Polyclonal (TR) Antibody - LS-C61136 - LSBio

CatalogID:	LS-C61136
Target:	Human IgG Fab'2
Host	Human IgG Fab'2 antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Texas Red (TR)
Immunogen Species:	Human IgG Fab'2 antibody was raised against Human
Immunogen:	Human IgG Fab'2 antibody was raised against sheep IgG F(ab')2 fragment.
Specificity:	anti-Rabbit Serum, Sheep IgG, Sheep IgG F(ab')2 and Sheep Serum. No reaction was observed against Sheep IgG F(c) or Human Serum Proteins
Reactivity:	Human
Purification:	Purified IgG
Presentation:	0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 10 mg/ml BSA (IgG, protease free), 0.05% sodium azide.
Recommended Storage:	Store at 4°C.
Usage Summary:	Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring lot-to-lot consistency.
Uses:	Immunofluorescence, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences