

Bovine IgG Rabbit anti-Bovine Polyclonal (Fab'2) (TR) Antibody - LS-C61262 - LSBio

CatalogID:	LS-C61262
Target:	Bovine IgG
Host	Bovine IgG antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Texas Red (TR)
Modifications:	Fab'2
Immunogen Species:	Bovine IgG antibody was raised against Bovine
Reactivity:	Bovine
Purification:	Purified IgG
Reconstitution:	deionized water. Possible additional volumes for resuspension: 500 µl
Uses:	Undetermined
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences