

Chicken IgY Rabbit anti-Chicken Polyclonal (Fab'2) (FITC) Antibody - LS-C61277 - LSBio

CatalogID:	LS-C61277
Target:	Chicken IgY
Host	Chicken IgY antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Modifications:	Fab'2
Immunogen Species:	Chicken IgY antibody was raised against Chicken
Antigen Type:	Purified protein
Immunogen:	Chicken IgY antibody was raised against chicken IgG whole molecule.
Specificity:	anti-Fluorescein, anti- Rabbit Serum, Chicken IgG and Chicken Serum. No reaction was observed against anti-Pepsin or anti-Rabbit IgG F(c)
Reactivity:	Chicken
Purification:	Immunoaffinity purified
Presentation:	0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 10 mg/ml chicken serum albumin (IgG, protease free), 0.01% sodium azide.
Recommended Storage:	Store at 4°C.
Usage Summary:	Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.
Uses:	Immunofluorescence (1:500 - 1:2500), Flow Cytometry (1:2000 - 1:10000) (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences