

Chicken IgY Rabbit anti-Chicken Polyclonal (Fab'2) (PE) Antibody - LS-C61281 - LSBio

<b>CatalogID:</b>	LS-C61281
<b>Target:</b>	Chicken IgY
<b>Host</b>	Chicken IgY antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Modifications:</b>	Fab'2
<b>Immunogen Species:</b>	Chicken IgY antibody was raised against Chicken
<b>Reactivity:</b>	Chicken
<b>Purification:</b>	Affinity purified
<b>Uses:</b>	Undetermined
<b>Size:</b>	500 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences