

Goat IgG Rabbit anti-Goat Polyclonal (Fab'2) (FITC) Antibody - LS-C61294 - LSBio

CatalogID:	LS-C61294
Target:	Goat IgG
Host	Goat IgG antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Modifications:	Fab'2
Immunogen Species:	Goat IgG antibody was raised against Goat
Antigen Type:	Purified protein
Immunogen:	Goat IgG antibody was raised against goat IgG whole molecule.
Specificity:	anti- Fluorescein, anti-Rabbit Serum, Goat IgG and Goat Serum. No reaction was observed against anti-Pepsin or anti- Rabbit IgG F(c). Limited reactivity will occur against Armenian Hamster IgG
Reactivity:	Goat
Purification:	Immunoaffinity purified
Reconstitution:	deionized water. Possible additional volumes for resuspension: 500 µl
Presentation:	0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 10 mg/ml goat serum albumin (IgG, protease free), 0.01% sodium azide.
Recommended Storage:	Store at 4°C.
Usage Summary:	Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.
Uses:	Immunofluorescence (1:500 - 1:2500), Flow Cytometry (1:2000 - 1:10000) (Optimal dilution to be determined by the researcher)
Size:	500 µg
Concentration:	1 mg/ml (after resuspension)
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences