

Pig IgG Rabbit anti-Pig Polyclonal (Fab'2) (TR) Antibody - LS-C61324 - LSBio	
<b>CatalogID:</b>	LS-C61324
<b>Target:</b>	Pig IgG
<b>Host</b>	Pig IgG antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Texas Red (TR)
<b>Modifications:</b>	Fab'2
<b>Immunogen Species:</b>	Pig IgG antibody was raised against Pig
<b>Reactivity:</b>	Pig
<b>Purification:</b>	Purified IgG
<b>Reconstitution:</b>	deionized water. Possible additional volumes for resuspension: 500 µl
<b>Uses:</b>	Undetermined
<b>Size:</b>	500 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	