

Syrian Hamster IgG Rabbit anti-Syrian hamster Polyclonal (Fab'2) (TR) Antibody - LS-C61333 - LSBio

CatalogID:	LS-C61333
Target:	Syrian Hamster IgG
Host	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Texas Red (TR)
Modifications:	Fab'2
Immunogen Species:	Syrian hamster
Antigen Type:	Purified protein
Immunogen:	Golden Syrian Hamster IgG whole molecule.
Specificity:	anti-Phycoerythrin, anti-Rabbit Serum, Golden Syrian Hamster IgG and Golden Syrian Hamster Serum. No reaction was observed against anti-Pepsin or anti-Rabbit IgG F(c)
Reactivity:	Syrian hamster, Hamster
Purification:	Purified IgG
Presentation:	0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 10 mg/ml BSA (IgG, protease free), 0.01% sodium azide.
Recommended Storage:	Store at 4°C. Do not freeze.
Usage Summary:	Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.
Uses:	Immunofluorescence (1:100 - 1:250), Flow Cytometry (1:100 - 1:250) (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences