

Horse IgG Rabbit anti-Horse Polyclonal (Fab'2) (PE) Antibody - LS-C61340 - LSBio

CatalogID:	LS-C61340
Target:	Horse IgG
Host	Horse IgG antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Modifications:	Fab'2
Immunogen Species:	Horse IgG antibody was raised against Horse
Reactivity:	Horse
Purification:	Affinity purified
Reconstitution:	Reconstitute with 1 ml of Deionized water.
Presentation:	Lyophilized in 0.02M Potassium Phosphate, 0.15M Sodium Chloride, pH 7.2
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Uses:	Undetermined
Size:	500 µg
Concentration:	0.5 mg/ml (after resuspension)
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences