

**Human IgG Donkey anti-Human Polyclonal (Fab'2) (TR) Antibody - LS-C61404 - LSBio**

<b>CatalogID:</b>	LS-C61404
<b>Target:</b>	Human IgG
<b>Host</b>	Human IgG antibody was produced in Donkey
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Texas Red (TR)
<b>Modifications:</b>	Fab'2
<b>Immunogen Species:</b>	Human IgG antibody was raised against Human
<b>Antigen Type:</b>	Purified protein
<b>Immunogen:</b>	Human IgG antibody was raised against human IgG whole molecule.
<b>Specificity:</b>	anti-Donkey Serum, Human IgG and Human Serum. No reaction was observed against anti-Pepsin, anti-Donkey IgG F(c) or Bovine, Chicken, Goat, Guinea Pig, Hamster, Horse, Mouse, Rabbit, Rat and Sheep Serum Proteins
<b>Reactivity:</b>	Human
<b>Purification:</b>	Immunoaffinity purified
<b>Reconstitution:</b>	deionized water. Possible additional volumes for resuspension: 500 µl
<b>Presentation:</b>	0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 10 mg/ml BSA, 0.01% sodium azide.
<b>Recommended Storage:</b>	Store at 4°C.
<b>Usage Summary:</b>	Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.
<b>Uses:</b>	Immunofluorescence (1:500 - 1:2500), Flow Cytometry (1:2000 - 1:10000) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	500 µg
<b>Concentration:</b>	1 mg/ml (after resuspension)
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences