

Rabbit IgG Goat anti-Rabbit Polyclonal (Fab'2) (TR) Antibody - LS-C61472 - LSBio

CatalogID:	LS-C61472
Target:	Rabbit IgG
Host	Rabbit IgG antibody was produced in Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Texas Red (TR)
Modifications:	Fab'2
Immunogen Species:	Rabbit IgG antibody was raised against Rabbit
Reactivity:	Rabbit
Purification:	Purified IgG
Reconstitution:	deionized water. Possible additional volumes for resuspension: 500 µl
Uses:	Undetermined
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences