

Rat IgG Fc Goat anti-Rat Polyclonal (Fab'2) (PE) Antibody - LS-C61516 - LSBio

CatalogID:	LS-C61516
Target:	Rat IgG Fc
Host	Rat IgG Fc antibody was produced in Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Modifications:	Fab'2
Immunogen Species:	Rat IgG Fc antibody was raised against Rat
Immunogen:	Rat IgG Fc antibody was raised against rat IgG Fc fragment.
Specificity:	anti-Phycoerythrin, anti-Goat Serum, Rat IgG, Rat IgG F(c) and Rat Serum. No reaction was observed against anti-Pepsin, anti-Goat IgG F(c), Rat IgG F(ab') ₂ or Bovine, Horse and Human Serum Proteins
Reactivity:	Rat
Purification:	Immunoaffinity purified
Presentation:	0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 10 mg/ml BSA (IgG, protease free), 0.01% sodium azide.
Recommended Storage:	Store at 4°C. Do not freeze.
Usage Summary:	Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.
Uses:	Immunofluorescence (1:100 - 1:250), Flow Cytometry (1:100 - 1:250) (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 8/29/2014

© 2014 LifeSpan BioSciences