

Rat IgG Goat anti-Rat Polyclonal (Fab'2) (TR) Antibody - LS-C61521 - LSBio

CatalogID:	LS-C61521
Target:	Rat IgG
Host	Rat IgG antibody was produced in Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Texas Red (TR)
Modifications:	Fab'2
Immunogen Species:	Rat IgG antibody was raised against Rat
Antigen Type:	Purified protein
Immunogen:	Rat IgG antibody was raised against rat IgG whole molecule.
Specificity:	anti-Goat Serum, Rat IgG and Rat Serum. No reaction was observed against anti-Pepsin, anti-Goat IgG F(c) or Bovine, Horse, Human, Mouse, Rabbit and Sheep Serum Proteins
Reactivity:	Rat
Purification:	Immunoaffinity purified
Presentation:	0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 10 mg/ml BSA (IgG, protease free), 0.01% sodium azide.
Recommended Storage:	Store at 4°C.
Usage Summary:	Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.
Uses:	Immunofluorescence, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences