

Rat IgG Rabbit anti-Rat Polyclonal (Fab'2) (TR) Antibody - LS-C61534 - LSBio

| | |
|-----------------------------|---|
| CatalogID: | LS-C61534 |
| Target: | Rat IgG |
| Host | Rat IgG antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Texas Red (TR) |
| Modifications: | Fab'2 |
| Immunogen Species: | Rat IgG antibody was raised against Rat |
| Antigen Type: | Purified protein |
| Immunogen: | Rat IgG antibody was raised against rat IgG whole molecule. |
| Specificity: | anti-Rabbit Serum, Rat IgG and Rat Serum. No reaction was observed against anti-Pepsin, anti-Rabbit IgG F(c) or Bovine, Horse, Human, Mouse, Rabbit and Sheep Serum Proteins |
| Reactivity: | Rat, Human |
| Purification: | Immunoaffinity purified |
| Presentation: | 0.01 M sodium phosphate, 0.15 M sodium chloride, pH 7.2, 10 mg/ml BSA (IgG, protease free), 0.01% sodium azide. |
| Recommended Storage: | Store at 4°C. |
| Usage Summary: | Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity. |
| Uses: | Immunofluorescence, Flow Cytometry (Optimal dilution to be determined by the researcher) |
| Size: | 500 µg |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences