

Rat IgG Donkey anti-Rat Polyclonal (Fab'2) (FITC) Antibody - LS-C61536 - LSBio

CatalogID:	LS-C61536
Target:	Rat IgG
Host	Rat IgG antibody was produced in Donkey
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Modifications:	Fab'2
Immunogen Species:	Rat IgG antibody was raised against Rat
Reactivity:	Rat
Purification:	Affinity purified
Uses:	Undetermined
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences