

Sheep IgG Donkey anti-Sheep Polyclonal (Fab'2) (TR) Antibody - LS-C61558 - LSBio

CatalogID:	LS-C61558
Target:	Sheep IgG
Host	Sheep IgG antibody was produced in Donkey
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Texas Red (TR)
Modifications:	Fab'2
Immunogen Species:	Sheep IgG antibody was raised against Sheep
Antigen Type:	Purified protein
Immunogen:	Sheep IgG antibody was raised against sheep IgG whole molecule.
Specificity:	anti- Donkey Serum, Sheep IgG and Sheep Serum. No reaction was observed against anti-Pepsin, anti-Donkey IgG F(c) or Chicken, Guinea Pig, Hamster, Horse, Human, Mouse, Rabbit and Rat Serum Proteins. Specificity was confirmed by ELISA at less than 1 % cross reactivity
Reactivity:	Sheep
Purification:	Immunoaffinity purified
Presentation:	0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 10 mg/ml BSA (IgG, protease free), 0.01% sodium azide.
Recommended Storage:	Store at 4°C.
Usage Summary:	Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.
Uses:	Immunofluorescence (1:500 - 1:2500), Flow Cytometry (1:2000 - 1:10000) (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences