

CD86 Rat anti-Mouse Monoclonal (PE) (RMMP-2) Antibody - LS-C62539 - LSBio

CatalogID:	LS-C62539
Target:	CD86 molecule
Synonyms:	CD86 Antibody, Activation B7-2 antigen Antibody, B7-2 Antibody, B7.2 Antibody, B70 Antibody, CD86 molecule Antibody, CTLA-4 counter-receptor B7.2 Antibody, CD86 antigen Antibody, FUN-1 Antibody, LAB72 Antibody, BU63 Antibody, CD28LG2 Antibody, Ctla-4 counter-receptors b7.2 Antibody
Host	CD86 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	RMMP-2
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD86 antibody was raised against Mouse
Specificity:	This antibody reacts with CD86 (B7-2) which is expressed on resting monocytes, dendritic cells, macrophages and at low levels on B lymphocytes ¹ . CD86 plays a role in T cell activation through interactions with CD28 and CTLA-4. Expression of CD86 can be en.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, EIA grade BSA.
Uses:	Flow Cytometry (1 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	50 µg or 300 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences