

**CD45R Rat anti-Mouse Monoclonal (FITC) (RA3-6B2) Antibody - LS-C62544 - LSBio**

<b>CatalogID:</b>	LS-C62544
<b>Target:</b>	CD45R
<b>Host</b>	CD45R antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	RA3-6B2
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CD45R antibody was raised against Mouse
<b>Immunogen:</b>	CD45R antibody was raised against mouse pre-B tumor cells (RAW112).
<b>Specificity:</b>	Anti-mouse CD45R, B220 (Ly 5) monoclonal antibody reacts with a form of the CD45 antigen found on B cells and lytically active subsets of NK cells and non - MHC restricted CTL's
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, 0.02% sodium azide, EIA grade BSA.
<b>Uses:</b>	Flow Cytometry (0.2 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg or 300 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences