

CD45R Rat anti-Mouse Monoclonal (PE) (RA3-6B2) Antibody - LS-C62545 - LSBio

CatalogID:	LS-C62545
Target:	CD45R
Host	CD45R antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	RA3-6B2
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD45R antibody was raised against Mouse
Immunogen:	CD45R antibody was raised against mouse pre-B tumor cells.
Specificity:	Anti-mouse CD45R, B220 (Ly 5) monoclonal antibody reacts with a form of the CD45 antigen found on B cells and lytically active subsets of NK cells and non - MHC restricted CTL's
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Uses:	Flow Cytometry (0.1 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	50 µg or 300 µg
Concentration:	0.1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences