

CD45RC Rat anti-Mouse Monoclonal (PE) (IBL-8) Antibody - LS-C62552 - LSBio

CatalogID:	LS-C62552
Target:	CD45RC
Host	CD45RC antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	IBL-8
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD45RC antibody was raised against Mouse
Antigen Type:	Cells
Immunogen:	CD45RC antibody was raised against mouse spleen cells.
Specificity:	CD45 (L-CA) is a transmembrane phosphotyrosine phosphatase expressed on leukocytes. Anti-mouse CD45RC is against the exon C-dependent RC isoform and reacts strongly with B cells, and less intensely with most CD8+ T cells. It does not recognize CD4+ T cell.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	Flow Cytometry (1 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	50 µg or 300 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences