

BCL3 / BCL-3 Hamster anti-Mouse Monoclonal (PE) (Ham150-3.5) Antibody - LS-C62565 - LSBio	
CatalogID:	LS-C62565
Target:	B-cell CLL/lymphoma 3 (BCL3)
Synonyms:	BCL3 Antibody, B-cell CLL/lymphoma 3 Antibody, B-cell leukemia/lymphoma 3 Antibody, B-cell lymphoma 3 protein Antibody, BCL-3 Antibody, BCL4 Antibody, D19S37 Antibody, Proto-oncogene BCL3 Antibody
Family / Subfamily:	Cancer
Host	BCL3 antibody was produced in Hamster
Clonality:	Monoclonal
Isotype:	IgG
Clone Name:	Ham150-3.5
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	BCL3 / BCL-3 antibody was raised against Mouse
Immunogen:	BCL3 / BCL-3 antibody was raised against recombinant Bcl-3-6-histidine tag.
Specificity:	Anti-mouse Bcl-3 monoclonal antibody reacts with both mouse and human Bcl-3 (B-cell Leukemia/lymphoma 3). The proto-oncoprotein Bcl-3 is a member of the I kappaB family and is present predominantly in the nucleus. This protein has been shown to function as.
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light. Store undiluted.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 µg or 300 µg
Concentration:	0.1 mg/ml
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	