

IL2RA / CD25 Rat anti-Mouse Monoclonal (FITC) (PC61.5.3) Antibody - LS-C62595 - LSBio	
CatalogID:	LS-C62595
Target:	interleukin 2 receptor, alpha (IL2RA)
Synonyms:	IL2RA Antibody, CD25 antigen Antibody, CD25 Antibody, IDDM10 Antibody, IL-2 receptor subunit alpha Antibody, IL2R Antibody, Interleukin 2 receptor, alpha Antibody, IL-2-RA Antibody, IL2-RA Antibody, TCGFR Antibody, IL-2R subunit alpha Antibody, p55 Antibody, TAC antigen Antibody
Family / Subfamily:	Interleukin
Host	IL2RA antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	PC61.5.3
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IL2RA / CD25 antibody was raised against Mouse
Antigen Type:	Cells
Immunogen:	IL2RA / CD25 antibody was raised against b6.1 CTL cell line.
Specificity:	Anti-mouse CD25 (IL-2R) monoclonal antibody reacts with the low affinity alpha chain of the interleukin-2 receptor antigen present on activated T and B cells in mice
Reactivity:	Mouse
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	Flow Cytometry (1 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	100 µg or 300 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	