

Ly76 / Ter119 Rat anti-Mouse Monoclonal (PE) (TER-119) Antibody - LS-C62599 - LSBio

CatalogID:	LS-C62599
Target:	Ly76 / Lymphocyte Antigen 76
Host	Ly76 / Lymphocyte Antigen 76 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	TER-119
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Ly76 / Ter119 antibody was raised against Mouse
Immunogen:	Ly76 / Ter119 antibody was raised against day 14 BALB/c fetal liver cells.
Specificity:	Anti-mouse erythroid cell (Ly-76) mAb is selectively reactive with both fetal and adult erythroid cells. This monoclonal antibody (clone: TER-119) is specific for cells at stages from early proerythroblast to mature erythrocytes. TER-119 is reported to re
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	Western blot, Flow Cytometry (0.05 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	50 µg or 300 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences