

**CD34 Rat anti-Mouse Monoclonal (FITC) (MEC14.7) Antibody - LS-C62602 - LSBio**

<b>CatalogID:</b>	LS-C62602
<b>Target:</b>	CD34 molecule
<b>Synonyms:</b>	CD34 Antibody, CD34 molecule Antibody, CD34 antigen Antibody
<b>Host</b>	CD34 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	MEC14.7
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CD34 antibody was raised against Mouse
<b>Specificity:</b>	Anti mouse CD34 reacts with Mouse CD34, a protein present on endothelial cells and hematopoietic progenitor cells. The antibody recognizes a neuraminidase sensitive epitope on endothelium in vivo, particularly on small vessels, neofomed capillaries and d.
<b>Reactivity:</b>	Mouse
<b>Uses:</b>	Undetermined
<b>Size:</b>	100 µg or 300 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences