

CD34 Rat anti-Mouse Monoclonal (PE) (MEC14.7) Antibody - LS-C62603 - LSBio

| | |
|---------------------------|--|
| CatalogID: | LS-C62603 |
| Target: | CD34 molecule |
| Synonyms: | CD34 Antibody, CD34 molecule Antibody, CD34 antigen Antibody |
| Host | CD34 antibody was produced in Rat |
| Clonality: | Monoclonal |
| Isotype: | IgG2a |
| Clone Name: | MEC14.7 |
| Conjugations: | Phycoerythrin (PE) |
| Immunogen Species: | CD34 antibody was raised against Mouse |
| Specificity: | PE anti-mouse CD34 monoclonal antibody. CD34 has been reported to be expressed on the endothelium of capillaries and, in this form, to function as a ligand for L selectin |
| Reactivity: | Mouse |
| Purification: | Purified |
| Presentation: | PBS, 0.1% sodium azide, EIA grade BSA. |
| Uses: | Flow Cytometry (1 µg/10E6 cells) (Optimal dilution to be determined by the researcher) |
| Size: | 50 µg or 300 µg |
| Concentration: | 0.1 mg/ml or 0.1 mg/ml |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences