

Sca-1 Rat anti-Mouse Monoclonal (PE) (CT-6A/6E) Antibody - LS-C62604 - LSBio

| | |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| CatalogID: | LS-C62604 |
| Target: | Sca-1 (Ly-6A.2/6E.1) |
| Host | Sca-1 (Ly-6A.2/6E.1) antibody was produced in Rat |
| Clonality: | Monoclonal |
| Isotype: | IgG2b |
| Clone Name: | CT-6A/6E |
| Conjugations: | Phycoerythrin (PE) |
| Immunogen Species: | Sca-1 antibody was raised against Mouse |
| Specificity: | Monoclonal antibody recognizes Sca-1 (Ly-6A.2/6E.1), a cell surface antigen used in the identification of hematopoietic stem cells. It is a member of the Ly-6 antigen family. The Thy-1lo, Lin- (lineage-negative, not expressing B220, Gr-1, Mac-1, CD4 or CD |
| Reactivity: | Mouse |
| Purification: | Purified |
| Presentation: | PBS, 0.02% sodium azide, EIA grade BSA. |
| Uses: | Flow Cytometry (2 µg/10E6 cells) (Optimal dilution to be determined by the researcher) |
| Size: | 50 µg or 300 µg |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences