

CD45RA Mouse anti-Rat Monoclonal (FITC) (MRC OX-33) Antibody - LS-C62629 - LSBio	
CatalogID:	LS-C62629
Target:	CD45RA
Host	CD45RA antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MRC OX-33
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD45RA antibody was raised against Rat
Immunogen:	CD45RA antibody was raised against rat purified spleen L-CA.
Specificity:	Rat Leukocyte Common Antigen has been shown to exist in different forms on different lymphoid cell types. Rat L-CA on thymocytes shows one main band at 180 kD, T cells show 4 bands at 180, 190, 200 and 220 kD. The antibody detects a subfraction of the 2.
Reactivity:	Rat
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	Flow Cytometry (0.5 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	100 µg or 500 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	