

CD2 Mouse anti-Rat Monoclonal (PE) (MRC OX-34) Antibody - LS-C62632 - LSBio

CatalogID:	LS-C62632
Target:	CD2 molecule
Synonyms:	CD2 Antibody, CD2 antigen Antibody, Erythrocyte receptor Antibody, LFA-3 receptor Antibody, LFA-2 Antibody, Lymphocyte-function antigen-2 Antibody, T11 Antibody, SRBC Antibody, CD2 molecule Antibody, Rosette receptor Antibody, T-cell surface antigen CD2 Antibody
Host	CD2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	MRC OX-34
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD2 antibody was raised against Rat
Immunogen:	CD2 antibody was raised against t blasts prepared in mixed lymphocyte reactions with purified rat T helper cells against irradiated spleen.
Specificity:	PE anti-rat CD2 (LFA-2) monoclonal antibody reacts with the 50 kD surface glycoprotein LFA-2, designated as CD2. LFA-2 is the receptor for LFA-3. This antibody labels all peripheral T cells and most thymocytes but does not label B cells or peritoneal mac.
Reactivity:	Rat
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	Flow Cytometry (0.1 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	50 µg or 200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences