

**Myeloid Cells Mouse anti-Rat Monoclonal (PE) (OX-82) Antibody - LS-C62643 - LSBio**

<b>CatalogID:</b>	LS-C62643
<b>Target:</b>	Myeloid Cells
<b>Host</b>	Myeloid Cells antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	OX-82
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	Myeloid Cells antibody was raised against Rat
<b>Specificity:</b>	Mouse anti-rat myeloid lineage monoclonal antibody recognizes a 35 kDa antigen found on myeloid cells and stromal elements from a variety of tissues in the adult rat. T
<b>Reactivity:</b>	Rat
<b>Uses:</b>	Western blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	50 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences