

IL2RA / CD25 Mouse anti-Rat Monoclonal (FITC) (MRC OX-39) Antibody - LS-C62648 - LSBio

CatalogID:	LS-C62648
Target:	interleukin 2 receptor, alpha (IL2RA)
Synonyms:	IL2RA Antibody, CD25 antigen Antibody, CD25 Antibody, IDDM10 Antibody, IL-2 receptor subunit alpha Antibody, IL2R Antibody, Interleukin 2 receptor, alpha Antibody, IL-2-RA Antibody, IL2-RA Antibody, TCGFR Antibody, IL-2R subunit alpha Antibody, p55 Antibody, TAC antigen Antibody
Family / Subfamily:	Interleukin
Host	IL2RA antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MRC OX-39
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IL2RA / CD25 antibody was raised against Rat
Immunogen:	IL2RA / CD25 antibody was raised against t blasts from a mixed lymphocyte reaction between purified CD4 positive T cells and irradiated spleens.
Specificity:	Anti-rat Interleukin-2 receptor monoclonal antibody recognizes the smaller (alpha subunit) 55kD chain of the IL-2 receptor found on activated rat T cells, thymic dendritic cells but not resting lymphocytes. The antibody binds to the rat interleukin -2 rece
Reactivity:	Rat
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	IHC - Paraffin, IHC - Frozen, Flow Cytometry (0.1 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	100 µg or 500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences