

SIRP / Signal-Regulatory Protein Mouse anti-Rat Monoclonal (PE) (MRC OX-41) Antibody - LS-C62652 - LSBio	
CatalogID:	LS-C62652
Target:	SIRP / Signal-Regulatory Protein
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	MRC OX-41
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Rat
Antigen Type:	Cells
Immunogen:	Rat peritoneal macrophages.
Specificity:	MRC OX-41 monoclonal antibody recognizes rat SIRP (Signal- Regulatory Protein), a surface protein of about 110 kDa. SIRP is selectively expressed by myeloid cells (macrophages, monocytes, granulocytes, dendritic cells) and neurons. The SIRP antigen is a t
Reactivity:	Rat
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	ICC, Immunofluorescence, Flow Cytometry (1 µg/10E6 cells), Radioimmunoassay (Optimal dilution to be determined by the researcher)
Size:	50 µg or 200 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	