

CD11b + c Mouse anti-Rat Monoclonal (PE) (MRC OX-42) Antibody - LS-C62655 - LSBio	
CatalogID:	LS-C62655
Target:	CD11b/CD11c
Host	CD11b/CD11c antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	MRC OX-42
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD11b + c antibody was raised against Rat
Antigen Type:	Cells
Immunogen:	CD11b + c antibody was raised against rat peritoneal macrophages.
Specificity:	Anti-rat CD1 1b/c monoclonal antibody recognizes most macrophages including resident peritoneal macrophages, Kupffer cells, but only approximately 35% of alveolar macrophages. The antibody also labels dendritic cells extensively, granulocytes and cells wi
Reactivity:	Rat
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Uses:	Flow Cytometry (0.2 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	50 µg or 200 µg
Concentration:	0.1 mg/ml
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	