

TCR Gamma+TCR Delta Mouse anti-Rat Monoclonal (FITC) (V65) Antibody - LS-C62669 - LSBio	
<b>CatalogID:</b>	LS-C62669
<b>Target:</b>	T Cell Receptor Gamma/Delta
<b>Host</b>	Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	V65
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	Rat
<b>Immunogen:</b>	Rat T Blasts.
<b>Specificity:</b>	Anti-rat gamma delta T cell receptor monoclonal antibody detects a T cell-specific heterodimeric 48 and 50 kD cell surface protein that is expressed on greater than 90% of CD3+ alpha beta TCR- rat peripheral T lymphocytes, and identifies a dense network o
<b>Reactivity:</b>	Rat
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 0.02% sodium azide, EIA grade BSA.
<b>Uses:</b>	Flow Cytometry (0.2 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg or 500 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	