

DPP4 / CD26 Mouse anti-Rat Monoclonal (PE) (OX-61) Antibody - LS-C62674 - LSBio

CatalogID:	LS-C62674
Target:	dipeptidyl-peptidase 4 (DPP4)
Synonyms:	DPP4 Antibody, ADCP-2 Antibody, ADABP Antibody, Dipeptidyl peptidase IV Antibody, Dipeptidyl-peptidase 4 Antibody, DPPIV Antibody, ADCP2 Antibody, T-cell activation antigen CD26 Antibody, CD26 Antibody, CD26 antigen Antibody, Dipeptidyl peptidase 4 Antibody, Dipeptidylpeptidase 4 Antibody, DPP IV Antibody, TP103 Antibody
Family / Subfamily:	Protease / Serine S9B
Host	DPP4 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	OX-61
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	DPP4 / CD26 antibody was raised against Rat
Immunogen:	DPP4 / CD26 antibody was raised against thoracic duct lymphocytes from a mesenteric lymphodectomized rat.
Specificity:	Anti-Rat CD26 monoclonal antibody (Clone: OX-61) recognizes and binds to the enzyme dipeptidyl peptidase IV (DPP-IV) present on rat thoracic duct lymphocytes and rat liver cells. The ectoenzyme DPP-IV is a bile canalicular cell surface molecule that has w.
Reactivity:	Rat
Purification:	Purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Recommended Storage:	Store at +4°C. DO NOT FREEZE.
Uses:	Flow Cytometry (0.5 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	50 µg or 200 µg
Concentration:	0.1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences